



Elasticated Cohesive Bandage

Weiyou® Elasticated Cohesive Bandage with excellent conformability and no sticky residue, it covers wounds and supports injured joints and muscles. As the manufacturer in China, Weiyou® has mature production lines and supply chains, and always supply best quality product, best favourable price and best service.

Product Description

Weiyou® Elasticated Cohesive Bandage is made of high quality porous fabric and non-woven material. The Elasticated Cohesive Bandage cohere only to itself, it does not adhere to other surfaces and provides good ventilation. The Elastic Cohesive Bandage is breathable, light, soft, less irritating to the skin, and can be quickly and easily applied to almost any body part, such as wrists, fingers, ankles, back, knees, feet and other body parts with your hands; it is also used for pets. The elastic cohesive bandage wrap is easy to tear by hand, no scissors needed.

Specification

Product name	Elastic Cohesive Bandage		
Material	Non-woven material		
Color	More than 10 different colors are available: white, green, blue, red, brown, etc		
Size	Size	Roll/box	Boxes/CTN
	1"×5yds	24	24
	2"×5yds	12	24
	3"×5yds	12	16
	4"×5yds	12	12
	4"×5yds	12	8
Logo	OEM customized logo is OK		

Certification	ISO / CE / FDA
Apply Occasion	Hospital, Clinic, Home

Features

Properties: Medical Adhesive & Suture Material

Size: 1"/2"/3"/4"/6"*5yds


Material: non-woven

Type: Medical Adhesive

Color: flesh,yellow,blue, green,red, black, etc







Strong Elastic



Self-Adhesive Style

No stick to hair, skin or clothing. Easy to tear.



Breathable Holes







SELF-ADHESIVE



HIGH ELASTIC





WEIYOU

Breathable

Waterproof

self-adhesive

Easy to tear

Latex free

HYPOALLERGENIC

Durable







